

XNAT Desktop: communicating non-DICOM data with XNAT

A Field Guide

Contents

- **XNAT Desktop: what, where, why**
- **Uploading files to XNAT**
- **File tagging techniques in XND**
- **Experimental features and future development**
- **Use of REST API by XND in communications with XNAT**

How do I upload non-DICOM data to XNAT?

- Manual
- Automatic script
- XNAT Desktop (XND)

XND: what

- **Cross platform**
- **Graphical user interface**
- **File-centric application**
- **Source code: Java+Eclipse Rich Client**
- **Beta version**

Attribute-value (tag) notion in XND

- Main XND objects: tags and files
- A tag is a <Name-Value> pair
- [tag1, tag2, ..] <-> individual file

Required tags to upload to XNAT

- XNAT container hierarchy:
Project << Subject << Experiment << [Scan, Assessment, Reconstruction]
- Modality: MR, CT, PET, US, ..
- Scan_Type: DICOM, TIFF, ...: required for XNAT image viewer
- (Optional): Quality, Date

XND Browsing

- Going through various file containers
- Manage/unmanage files
- Assigning tags to files within containers


Folder containers in XND: File System

- Managed root
- Special: <Incoming> folder
- Subfolders under managed root


Non-FS file containers

- **Virtual Folder (VF):** certain tags can be assigned to show as VF's in hierarchy
- **Collection (~XNAT resource):** all files within collection share the same tags.

Directory browsing



Tagged browsing



WAYS TO SET TAGS

1. Manual entry
2. Copy & paste
3. Excel spreadsheet
4. Automatic tag extraction rules

XNAT TAGS

Project	sbp	Subject	subj1
Experiment	mr1	Scan	1
Time		Date	
Description			
USER TAGS			
measure1	22	measure2	1.24

XND folder view

The screenshot displays the XNAT interface. On the left, a 'File view' pane shows a hierarchical folder structure. The 'data' folder is expanded, revealing sub-folders like 'qBOLD', 'Control_qBOLD_CCIR_EV005...', 'DICOM', 'map_t2.dcm', 'rqb_out', and 'Gamma_qBOLD' series. On the right, a table lists resources with columns for 'Resource Name', 'Project', 'Subject', 'Experiment', and 'Scan'. The selected resource is 'default_EV005_CCIR_02...', with a 'scan_9_catalog.xml' file icon next to it. Below the table, an 'All tags' section shows a list of metadata tags and their values.

Resource Name	Project	Subject	Experiment	Scan
default_EV005_CCIR_02...	default	EV005	CCIR_02...	9

Name	Value
Collection_ID	data_qBOLD_Control_qBOLD_CCIR_EV005_...
Date	20090417
Experiment	CCIR_0209_Mintun_Benzinger_TOF_MRA_...
Modality	MR
Project	default
Scan	9
SeriesDescription	gre_field_mapping_for_OEF

XND tag view

The screenshot displays the XNAT interface in 'Tag view' mode. On the left, a tree view shows the hierarchy: Project:default > Subject:EV005 > Experiment:CCIR_0209 > Scan:9. The right pane shows a table with the following data:

Resource Name	Project	Subject	Experiment	Scan
default_EV005...	default	EV005	CCIR_0209_Mintun...	9

Below the table, the 'All tags' section contains the following information:

Name	Value
Collection_ID	data_qBOLD_Control_qBOLD_CCIR_EV005_map_t2.d
Date	20090417
Experiment	CCIR_0209_Mintun_Benzinger_TOF_MRA__QBOLD_
Modality	MR
Project	default
Scan	9

Upload from XND to XNAT

- Use Tag View to upload
- Can upload files at any level
- Single/multiple container upload

(Upload demo)

Steps to tag files in XND

- Add managed root (rule of thumb: use lowest level directory possible)
- Manage some files in subfolders
- Select a container
- Apply tagging procedure to the container
- Repeat using tag copy/paste feature

Manual tagging example (see the demo)

- Manual: set 'Project' and 'quality'
- Manual: add Experiment and Scan values to one of subfolders
- Use clipboard copy/paste to quickly populate tags

Sources of information about local files

- File contents
- File names & extensions
- Directory path
- External sources (spreadsheets, xml, etc.)

Automatic tagging: rules

- Metadata source
- Default logic (built-in),
- Custom logic (user-defined in xml file using rule language)

Automatic tagging: DICOM rule

- **Source: DICOM tags**
- **Default logic: create collections based on matching tag hierarchies**
- **Custom logic: default (XNAT-specific) or user-supplied tag map**

DICOM rule tag matching

- (element, group) <-> XND tag
- Default value
- Assignment prioritization


```
<tag name="Subject" defaultValue="subj_undefined">  
  <DICOMTaggroup="0010" element="0010" alias="PatientName" priority="1" />  
  <DICOMTaggroup="0010" element="0020" alias="PatientID" priority="2" />  
</tag>
```

(DICOM rule demo)

Automatic tagging: naming rule


- **Metadata source: folder structure/file name pattern**
- **Default logic: recursive processing of selected folder hierarchy**
- **Custom logic: match parts of file/folder names with regular expressions**

Naming rule: modeling folder structure




```
<folderID="root_fold" pattern=".*" treeRoot="1">
  <child>AnyFolder </child>
  <child>label </child>
</folder>
<folderID="AnyFolder" pattern=".*">
  <tag name="FolderName" pattern=".*" recursive="pattern" />
  <child>AnyFolder </child>
</folder>
```

Naming rule: capturing folder name patterns




```
<folderID="OAS1" pattern="OAS1.*" treeRoot="1">  
  <tag name="Subject" pattern="(.{9}).*" recursive="fixed" />  
</folder>  
<folderID="FSL_SEG">
```


Naming rule: file name patterns



Folder ID =
Folder Name



```
<folderID="RAW">  
  <child>SCAN </child>  
</folder>  
<fileID="SCAN" pattern=".*">  
  <tag name="Scan" pattern=".*_mpr-(.)_anon.*" recursive="pattern" />  
</file>
```

(Naming rule demo)

Automatic tagging, tag pattern rule

- **Source:** existing tags in XND
- **Default logic:** use specified tag to derive a modified value
- **Custom logic:** user-defined tag map (regexp based)

Tag pattern: formatting subject name example

```
<tag name="Subject">  
  <substring pattern=".{4}(.{3})*" value="0{1}"/>  
</tag>  
<tag name="Scan_Type">  
  <replace match="[\x3c\x5b\x3e\x5d]" with="" />  
  <replace match="[" with="_" />  
  <replace match="[-]" with="_" />  
  <replace match="[\x26]" with="_" />  
</tag>
```

file1 [Subject="wcaf01", Scan_Type="pre-dti.2mm&low^contrast"]



file1 [Subject="001", Scan_Type="pre_dti_.2mm_low_contrast"]

Macro: a tool for one-click data markup

- Macro is a sequence of commands of three following types:
 - Manage
 - Apply rule
 - Set a tag
- Construct macros for re-occurring data patterns for quick tag & upload.

(Macro demo)

Conclusions

- Use XNAT Desktop for uploading multi-format data to XNAT
- Use supplied tagging rules or create your own
- Tagging procedure can, with some customization, be made automatic in most cases

What's next? (Features in development. Early previews available)

- Data import wizard
- DICOM networking capabilities: import directly from PACS
- Research PACS workstation interface: extended DICOM image markup and viewing

XND: building REST queries

- XNAT:

Project << Subject << Experiment << [Scan,
Assessment, Reconstruction


- XND:

<Project="default"><Subject="EV005"><Experiment="CCIR00209.."><Scan="9">

Tag relationships in XND

```
<ontology id="xnat" descr="Default ontology with XNAT hierarchical tags.">
  <tag name="Project" treeRoot="1" tableView="1" anValue="1" type="ontology">
    <child>Subject</child>
    <child>Resource</child>
  </tag>
  <tag name="Subject" tableView="1" anValue="1" type="ontology" context="1">
    <child>Experiment</child>
    <child>Resource</child>
  </tag>
  <tag name="Experiment" tableView="1" anValue="1" type="ontology">
    <child>Scan</child>
    <child>Assessor</child>
    <child>Reconstruction</child>
    <child>Resource</child>
  </tag>
```

XND <-> XNAT communication



XND upload strategy

- Determine level at which upload should be performed.
- Execute GET query at each hierarchy level. Find lowest non-existent level.
- Re-create hierarchy using PUT methods that create XNAT objects
- Upload under created hierarchy.
- Repeat for multiple items

REST upload example

GET (success): [http://central.xnat.org/REST/experiments?](http://central.xnat.org/REST/experiments?xsiType=xnat:subjectAssessorData&project=default&format=xml)

[xsiType=xnat:subjectAssessorData&project=default&format=xml](http://central.xnat.org/REST/experiments?xsiType=xnat:subjectAssessorData&project=default&format=xml)

GET (error code 404):

http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/scans/9/?format=xml

GET (error code 404): http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/?format=xml

GET (success): <http://central.xnat.org/REST/projects/default/subjects/EV005/?format=xml>

PUT (success): [http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/?](http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/?date=20090417&label=CCIR_0209&xsiType=xnat:MRSessionData&format=xml)

[date=20090417&label=CCIR_0209&xsiType=xnat:MRSessionData&format=xml](http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/?date=20090417&label=CCIR_0209&xsiType=xnat:MRSessionData&format=xml)

PUT (success): [http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/scans/9/?](http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/scans/9/?ID=9&series_description=gre_field_mapping_for_OEF&format=xml)

[ID=9&series_description=gre_field_mapping_for_OEF&format=xml](http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/scans/9/?ID=9&series_description=gre_field_mapping_for_OEF&format=xml)

PUT (success):

[http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/scans/9/resources/DICOM/?](http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/scans/9/resources/DICOM/?content=RAW&format=DICOM&)

[content=RAW&format=DICOM&](http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/scans/9/resources/DICOM/?content=RAW&format=DICOM&)

PUT (success):

[http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/scans/9/resources/DICOM/files/fil](http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/scans/9/resources/DICOM/files/files.zip?inbody=true&extract=true&format=DICOM&content=RAW&)

[es.zip?inbody=true&extract=true&format=DICOM&content=RAW&](http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/scans/9/resources/DICOM/files/files.zip?inbody=true&extract=true&format=DICOM&content=RAW&)

PUT (success): [http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/?](http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/?pullDataFromHeaders=true&format=xml)

[pullDataFromHeaders=true&format=xml](http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/?pullDataFromHeaders=true&format=xml)

PUT (success): [http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/?](http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/?triggerPipelines=true&format=xml)

[triggerPipelines=true&format=xml](http://central.xnat.org/REST/projects/default/subjects/EV005/experiments/CCIR_0209/?triggerPipelines=true&format=xml)

Resources

- XND manual:
<http://www.xnat.org/xnd/>
- XND download:
<http://nrg.wustl.edu/xnd/download/>
- Java regular expression specification:
<http://java.sun.com/j2se/1.4.2/docs/api/java/util/regex/Pattern.html>
- Regular expression test applet:
<http://www.cis.upenn.edu/~matuszek/General/RegexTester/regex-tester.html>