

DICOM services for XNAT

Using the DICOM standard to streamline
data storage and retrieval


Plan


In this talk, we will:

- See how users can get DICOM data into and out of XNAT
- Identify the components XNAT uses for DICOM integration
- Learn how to configure and troubleshoot XNAT's DICOM services

Introduction

- XNAT DicomServer simplifies data import
- DicomBrowser offers GUI and command-line tools for finding, deidentifying, and sending DICOM data
- XNAT Gateway provides DICOM query and retrieve services





Study Description	Study Date and Time	Series	Instances
COR00100 COR-00101 MARKS	6/20/10 12:38:15 PM	22	1346




TAZINTOSH.COM

Edward Pyra


Components

- DicomServer uses the XNAT web application as a store for session metadata
- Imaging data are stored as files on disk, not in the XNAT database




Components

- **DicomServer:**
 - is a DICOM C-STORE receiver (SCP)
 - receives DICOM data from a network and saves files to disk
 - translates DICOM metadata to XNAT metadata format
- XNAT can also save files and extract metadata




Components: DicomServer

Any DICOM C-STORE client can send data, including scanners...



Components: DicomServer


...also dedicated workstations, or viewers like OsiriX or DicomBrowser...



Components: DicomServer

...or the XNAT upload applet, which:

- guides user to label the XNAT session
- uploads DICOM or ECAT data



Uploading data

Upload client needs to know some connection parameters:

- server hostname
- port number
- DICOM AE title

Launch Uploader

Option 3: DICOM Browser (via DICOM Server)

Use a custom built DICOM application which supports easy anonymization, and can perform a DICOM SEND to a specified DICOM Server.

DICOM SERVER Specifications

- Host Name: cnda.wustl.edu
- Port: 104
- AE Title: CNDA

Start DICOM Browser

powered by
XNAT

Uploading data

Which XNAT project receives the data can be specified in:

- Patient Comments
- Study Comments
- Study Description
- Accession Number
- or some other (configurable) DICOM attribute

Uploading data

Received data (usually) goes into XNAT prearchive

- Prearchive is a holding/review area
- Projects can be configured to skip prearchive step (“autoarchiving”)
- Autoarchiving can be specified per session in DICOM Patient Comments or Study Comments as `AA:true` or `AA:false`

XNAT upload applet

- Provides a GUI for finding, labeling, and uploading DICOM data (ECAT too)
- Uploader uses DICOM C-STORE to move DICOM data
- Also uses second channel (HTTP) to carry XNAT session metadata


Data model conversion

DicomServer builds an XNAT session from each DICOM study it receives

- Data are reorganized into files in XNAT directory structure
- Session directory is placed in appropriate location (XNAT prearchive or archive)
- DICOM metadata are translated into XNAT session and scan information

Uploading, in detail



- DICOM C-STORE to transfer data
- No explicit end-of-study signal; DicomServer uses idle timeout to decide when a session is complete
- DicomServer translates DICOM metadata into XNAT format



Metadata conversion

XNAT and DICOM have different metadata models

- DICOM study ~ XNAT session
- DICOM series ~ XNAT scan



MR Session: Sample_ID

Details Projects

Accession #	Tyto_E00001	Subject:	Sample_Patient
Date Added	2010-05-13 12:27:24.0 (karchie)	Gender:	
Date:	2006-12-14	Handedness:	
Time:	09:12:06	Age:	--
Scanner:	MEDPC SIEMENS TrioTim		
Acquisition Site:	Hospital		

Notes:

Scans

Scan	Type	Usability	Files
4	t1_mpr_1mm_p2_pos50	usable	DICOM (176 files, 32.21 Mb) SNAPSHOTS (2 files, 343 Kb)
5	t1_mpr_1mm_p2_pos50	usable	DICOM (176 files, 32.11 Mb) SNAPSHOTS (2 files, 343 Kb)
6	t2_spc_1mm_p2	usable	DICOM (176 files, 32.38 Mb) SNAPSHOTS (2 files, 277 Kb)

Metadata conversion

Each XNAT metadata field is derived from one or more DICOM fields

```
final class ImageScanAttributes {  
    private ImageScanAttributes() {} // no instantiation  
  
    static public ReadableAttrDefSet<Integer,String> get() { return s; }  
  
    static final private AttrDefSet s = new AttrDefSet();  
  
    static {  
        s.add("ID"); // handled by session builder  
        s.add("UID", Tag.SeriesInstanceUID);  
  
        s.add("series_description", Tag.SeriesDescription);  
        s.add("scanner", Tag.StationName);  
        s.add("scanner/manufacture", Tag.Manufacturer);  
        s.add("scanner/model", Tag.ManufacturerModelName);  
    }  
}
```

Metadata conversion

Some conversions are nontrivial and require their own Java classes

```
s.add(new XnatAttrDef.Time("parameters/scanTime", Tag.SeriesTime));
s.add(new VoxelResAttribute("parameters/voxelRes"));
s.add(new OrientationAttribute("parameters/orientation"));
s.add("coil", Tag.ReceiveCoilName);
s.add("fieldStrength", Tag.MagneticFieldStrength);
s.add(new XnatAttrDef.Real("parameters/tr", Tag.RepetitionTime));
s.add(new MREchoTimeAttribute());
s.add(new XnatAttrDef.OptionalWrapper(new XnatAttrDef.Real("parameters/ti", Tag.InversionTime)));
s.add(new XnatAttrDef.Int("parameters/flip", Tag.FlipAngle));
s.add("parameters/sequence", Tag.SequenceName);
s.add("parameters/imageType", Tag.ImageType);
s.add("parameters/scanSequence", Tag.ScanningSequence);
s.add("parameters/seqVariant", Tag.SequenceVariant);
s.add("parameters/scanOptions", Tag.ScanOptions);
s.add("parameters/acqType", Tag.MRAcquisitionType);
s.add(new ImageFOVAttribute("parameters/fov"));
```

Metadata modification

Upload applet and DicomBrowser have a scripting language for specifying changes to DICOM metadata (a.k.a. “anon scripts”)

```
(0010,0010) := "Subject_346"  
-(0010,0030)          // Patient's Birth Date  
(0008,0080) = "CCIR 3T" : (0008,0080) := "WUSM CCIR"  
(0020,000d) := new UID  
(0010,0020) := replace[(0010,0020), " ", "_"]
```

Metadata modification

Each XNAT project may be configured with a script that is applied to DICOM data by the upload applet:

```
http://my.xnat.org/REST/projects/MyProject/  
resources/UPLOAD_CONFIG/files/dicom.das
```


Configuration: basic

```
<?xml version="1.0" encoding="UTF-8"?>
<!DOCTYPE properties SYSTEM "http://java.sun.com/dtd/properties.dtd>
<properties>
<comment>Simple DicomServer configuration</comment>

<entry key="arcspec">/opt/xnat/cache/archive_specification.xml</entry>
<entry key="xnat_url">http://my.xnat.org:8080/xnat</entry>
<entry key="user">admin</entry>
<entry key="password">adminpw</entry>
<entry key="http_port">8180</entry>

<entry key="log4j.rootLogger">WARN,R</entry>
    ...
</properties>
```

Configuration:port mapping

Can map port 104 to 8104 (via iptables) without changing DicomServer config

DICOM host

This is the host where the C-Store is accessible. It should match the value you set when you configured DICOM Server as 'public_host'. This is not a URL and should not include a scheme (http or https). It should just be the host name (the name of the machine on which it is running).

DICOM port

This is the port where the C-Store is accessible. It should match the value you set when you configured DICOM Server as 'public_port'. DICOM Server's default value is 104.


DICOM AE Title

This is the port where the C-Store is accessible. It should match the value you set when you configured DICOM Server as 'ae_title'. DICOM Server's default value is XNAT.

DICOM Server HTTP URL

This is the url used by the applet to connect to DICOM Server's internal HTTP server. It should match the value you set when you configured DICOM Server as 'sessions_url'.

Configuration:reverse proxy



Configuration:reverse proxy

```
<entry key="arcspec">/opt/xnat/cache/archive_specification.xml</entry>
<entry key="xnat_url">http://my-hidden-host.local:8080/xnat</entry>
<entry key="user">admin</entry>
<entry key="password">adminpw</entry>

<entry key="http_port">8180</entry>
<entry key="xnat_public_url">https://my.xnat.org</entry>
<entry key="scp_public_host">my.xnat.org</entry>
<entry key="scp_public_port">104</entry>
<entry key="sessions_url">https://my.xnat.org:8443/session</entry>

<entry key="log4j.rootLogger">WARN,R</entry>

...
```


Troubleshooting

- Three places to find diagnostic info:
- DicomServer log (and received file log)
- upload applet Java console
- session management web service



Troubleshooting FAQs

Where did my upload go?

- If DicomServer can't determine to which project a session should go, it goes to the Unassigned prearchive
- Only admin accounts can view the Unassigned prearchive
- If the session isn't in Unassigned, check main DicomServer log and received log

Troubleshooting

Why didn't my session autoarchive?

Common reasons:

- The session is already archived in that project (Study Instance UID)
- The session label is already used in that project
- The project isn't set to autoarchive

Troubleshooting

Why did my study get split?

- DICOM C-STORE protocol doesn't carry any information about study size
- DicomServer has to guess when the sender has finished with a study
- Once DicomServer starts autoarchiving a study, newly received files go into the prearchive

Summary

- XNAT uses DicomServer to receive DICOM data
- Users shouldn't (and usually don't) need to know about DicomServer



Links

DicomServer:

<http://www.xnat.org/DicomServer>

DicomBrowser:

<http://nrg.wustl.edu/projects/DICOM/DicomBrowser.jsp>

<http://nrg.wustl.edu/projects/DICOM/DicomBrowser-cli.html>

XNAT Gateway:

<http://www.xnat.org/XNAT-DICOM+Gateway>

Questions?

Photo credits:
Edouard Puginier
Scott Moynihan
Tony Case
Daniel Rothamel
Steve Wall