

Creating and Managing Projects

Data stored within XNAT is organized according to projects. To manage data within your XNAT you will need to create projects which your data can be assigned to. In some situations, these projects correspond directly to grant approved studies. In other situations, projects can correspond to individual collection sites for a multi-site study. You are free to use any project organization scheme you choose. However, remember that user permissions for access to data are governed by the user's access to particular projects. So, if you only have one project in your XNAT, you will have limited options for governing access to that data.

All data stored in XNAT must be associated with a project. This association is the basis of the XNAT security model. Users are given access to data which belongs to a particular project. In order to upload any data to XNAT, you must first create a project.

Creating a Project

Go to New -> Project in the top navigation

New Project

Step 1: Enter project details

Project Title

REQUIRED: Enter the full name of your project here. This will show up on project listings.

Running Title

REQUIRED: Create a simple abbreviation of your project name, using 24 characters or less. Spaces are allowed. This will be commonly used in menus and UI elements.

Project ID

REQUIRED: Create a one word project identifier. This is used in the database and **cannot be changed**.

Project Description

Optional: Provide a description of your project. This is for reference only and is not searchable.

Keywords

Optional: Enter searchable keywords. Each word, separated by a space, can be used independently as a search string.

Alias(es)

Optional: Enter alternate aliases (for example: charge codes) that this project can be identified by.

Investigator(s) **Primary Investigator:** [Edit Selected](#)

Other Investigators:

Optional: List investigators associated with this project. This is for reference only and **does not provide access** to this project for the listed investigators.

Required Fields:

Project Title	The full title of your project, which will appear in the site-wide project listing.
Running Title	This is an abbreviated project name, which will appear in search listings and data tables.
Project ID	This ID is used within the XNAT database and in all API calls related to your project and its data. Once set, it cannot be changed.

Optional Fields:

Description	Add as much detail as you like describing your project. This will show up in the project page, and in the site-wide project listing.
Keywords	Keywords can be searched for in the primary search field, and can help you find project data in an XNAT system with many projects.
Aliases	These are optional project identifiers. A common reason to use them would be for billing codes.

Setting Project Investigators



Assigning an investigator to a project does not grant that person access to the project. A user account is required to gain project access.

Project Investigators can be added to projects as reference values only. These do not need to correspond to user accounts in your XNAT webapp. A common reason to use this would be for identification of primary investigators or researchers associated with a study, or with multiple studies.

You can register investigators in your XNAT by clicking the "Create Investigator" button in the Add Project form.

You can associate any registered investigators in your XNAT with your project. You can also assign a single investigator as a primary investigator.

The screenshot shows a 'Create Investigator' dialog box with the following fields and options:

- Title: [Text input field]
- First Name: [Text input field]
- Last Name: [Text input field]
- Institution: [Text input field]
- Department: [Text input field]
- Email: [Text input field]
- Phone: [Text input field]
- Set as Primary
- Buttons: Cancel, Submit

Setting Project Accessibility

This setting allows you to manage who can view project data, or even have awareness of the existence of your project in XNAT. There are three possible accessibility settings.

Public	All XNAT users will be able to view or download your project and its data. Owners of other XNAT projects can share your data into their studies. If your XNAT security settings allow guest access, then any visitor to the web application will be able to view and download your data.
---------------	--

Protected	All XNAT users will be able to see your project title and description in site-wide listings, but only users assigned by project owners will be able to access its data.
Private	Only XNAT users assigned by project owners will be able to see the project in XNAT, or access its data.

Saving Project Data

After completing the project creation dialog, you will be taken to the Project report page. This page is the central place for dealing with your new project. From here you can manage user access, manage project protocol details, monitor project progress, and directly manage stored data. Take some time to familiarize yourself with the options the Project Report page offers. They are different depending on your role in the project (Project Owners can manage user access, Project Members can create and edit data, Project Collaborators can only read/download data).

Managing User Access

Unless you have set your project to "Public" in the Project Accessibility setting above, access to your project data will be controlled by users you have specified in the "Access" tab on the Project page.

Dicom to Nifti Conversion Project

Details **Access** Manage Pipelines

Show Deactivated Users

Username	First Name	Last Name	Email	Group	Remove
<input type="text" value="Filter by Username"/>	<input type="text" value="Filter by First N"/>	<input type="text" value="Filter by Last N"/>	<input type="text" value="Filter by Email"/>		
admin	Admin	Admin	administrator@xnat.org	Owners	✗
owner	Sample	User	xnatselenium@gmail.com	Owners	✗
member	Sample	User	xnatselenium@gmail.com	Members	✗
collaborator	Sample	User	xnatselenium@gmail.com	Collaborators	✗

Add/Invite User

Enter the username or email address of the user(s) to add. Separate multiple users by spaces or commas.

Add Users from List

Project Access Settings

Define Default Project Accessibility

Private Only you and study members will be able to access study data.

Protected All users will be able to see your study title and description, but only collaborators you approve will be able to analyze and download data.

Public All users will be able to access study data for analysis and download.

Project Invitations

No outstanding invitations to this project

You can invite existing users by clicking the "Show List" button, or you can invite new users by their email address by clicking the "Add User" button. In either case, you have to choose what access level to give them: Project Owners, Project Members or Project Collaborators.

See: [Understanding User Roles and Permissions](#).

Data Ownership and Sharing

Subjects and experiments are always owned by one (and only one) project. This project governs the ability of users to modify the data. However, it is common for a user with access to your project to share its data into their own project. When this occurs, your data cannot be modified in this second project, and any modification made to the data in your original project is also reflected in the second project.

See: [Understanding Data Sharing in XNAT's Security Structure](#).

Project Reporting

The progress of your project can be monitored from the Project Report. The report includes two important sections for monitoring progress: Details & Listings. The Listings section displays detailed information about the data you have stored in your project.

The Project Listings section is a central point for XNAT users. From it you can quickly view the data within a project and link to particular entries. By default, data tables for Subjects and MR Sessions will be preloaded. However, you can add or remove data tables for any data type that exists in your project. You can also create and display any project-specific saved search here.

See: [Using Stored Searches](#)