

Experiment Data REST XML Path Shortcuts

xnat:experimentData

These URIs require the use of the **xsiType** parameter if you are using the query string option. The following parameters are available at these URIs:

```
/data/projects/{ID}/subjects/{ID}/experiments  
/data/projects/{ID}/subjects/{ID}/experiments/{ID}
```

Experiment Data Parameters

Variable Name	XML Path
visit_id	xnat:experimentdata/visit_id
date	xnat:experimentdata/date
ID	xnat:experimentdata/ID
project	xnat:experimentdata/project
label	xnat:experimentdata/label
time	xnat:experimentdata/time
note	xnat:experimentdata/note
pi_firstname	xnat:experimentdata/investigator/firstname
pi_lastname	xnat:experimentdata/investigator/lastname
validation_method	xnat:experimentdata/validation/method
validation_status	xnat:experimentdata/validation/status
validation_date	xnat:experimentdata/validation/date
validation_notes	xnat:experimentdata/validation/notes
last_modified (READ-ONLY)	xnat:experimentdata/meta/last_modified
insert_date (READ-ONLY)	xnat:experimentdata/meta/insert_date
insert_user (READ-ONLY)	xnat:experimentdata/meta/insert_user

Subject Data Parameters

Variable Name	XML Path
subject_ID	xnat:subjectData/ID
subject_label	xnat:subjectData/label
subject_project	xnat:subjectData/project

Image Session Data Parameters

Variable Name	XML Path
scanner	xnat:imageSessionData/scanner
operator	xnat:imageSessionData/operator
dcmAccessionNumber	xnat:imageSessionData/dcmAccessionNumber
dcmPatientId	xnat:imageSessionData/dcmPatientId
dcmPatientName	xnat:imageSessionData/dcmPatientName
session_type	xnat:imageSessionData/session_type

modality	xnat:imageSessionData/modality
UID	xnat:imageSessionData/UID

MR Session Data Parameters

Variable Name	XML Path
coil	xnat:mrSessionData/coil
fieldStrength	xnat:mrSessionData/fieldStrength
marker	xnat:mrSessionData/marker
stabilization	xnat:mrSessionData/stabilization

PET Session Data Parameters

Variable Name	XML Path
studyType	xnat:petSessionData/studyType
patientID	xnat:petSessionData/patientID
patientName	xnat:petSessionData/patientName
stabilization	xnat:petSessionData/stabilization
scan_start_time	xnat:petSessionData/start_time/scan
injection_start_time	xnat:petSessionData/start_time/injection
tracer_name	xnat:petSessionData/tracer/name
tracer_startTime	xnat:petSessionData/tracer/startTime
tracer_dose	xnat:petSessionData/tracer/dose
tracer_sa	xnat:petSessionData/tracer/specificActivity
tracer_totalmass	xnat:petSessionData/tracer/totalMass
tracer_intermediate	xnat:petSessionData/tracer/intermediate
tracer_isotope	xnat:petSessionData/tracer/isotope
tracer_isotope	xnat:petSessionData/tracer/isotope/half-life
tracer_transmissions	xnat:petSessionData/tracer/transmissions
tracer_transmissions_start	xnat:petSessionData/tracer/transmissions/startTime